

ANGULATED DISSECTING FORCEPS 16 CM

Ref. IP04516



07/04/2026 01:25:48 – For information purposes only, can be subject to change

ANGULATED DISSECTING FORCEPS 16 CM

Ref. [IP04516](#)



Allows to carry out care and facilitate the gripping of small objects

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021843	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.026	Bare Unit / Weight (kg)	0.025

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 01:25:48 – For information purposes only, can be subject to change