

DISSECTING SET

Ref. IT01600



07/04/2026 01:19:34 – For information purposes only, can be subject to change

DISSECTING SET

Ref. IT01600



Complete dissection kit

Description

Consisting of:

- IP26410 Splinter forceps, 10 cm
- IP26014 Dissection forceps without teeth, 14 cm
- IS02414, Fluted probe, 14 cm
- ID07601, Dissecting needle with handle, 14 cm
- ID08601, Dissecting cutting needle with handle, 14 cm
- IC11414 Straight round-ended scissors, 14 cm
- IC09411 Straight iridectomy scissors, 11.5 cm
- ID05404, no.4 scalpel handle

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

DISSECTING SET

Ref. IT01600



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040031569	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	25	Bare Unit / Width (cm)	13.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	15	Bare Unit / Length (cm)	18
Packed unit / Weight	0.282	Bare Unit / Weight (kg)	0.282
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 01:19:34 – For information purposes only, can be subject to change