

**S/S KOCHER FORCEPS JAR 25 (DIAM)X H.
140 MM**

Ref. AI11615



07/04/2026 00:02:08 – For information purposes only, can be subject to change

**S/S KOCHER FORCEPS JAR 25 (DIAM)X H.
140 MM**



Ref. AI11615

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019307	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	9	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	17	Bare Unit / Length (cm)	14
Packed unit / Weight	0.126	Bare Unit / Weight (kg)	0.121
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 00:02:08 – For information purposes only, can be subject to change