

**S/S SERVING FORCEPS JAR 200 X 50 (DIAM)
MM**

Ref. AI11623



CE



07/04/2026 01:42:30 – For information purposes only, can be subject to change

**S/S SERVING FORCEPS JAR 200 X 50 (DIAM)
MM**



Ref. AI11623

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019314	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	10	Bare Unit / Width (cm)	5
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20
Packed unit / Weight	0.213	Bare Unit / Weight (kg)	0.213
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 01:42:30 – For information purposes only, can be subject to change