

# S/S BOWL WITH SPOUT DIAM 240 MM - 3LITERS

Ref. AI15240



07/04/2026 02:38:13 – For information purposes only, can be subject to change

# S/S BOWL WITH SPOUT DIAM 240 MM - 3LITERS

Ref. AI15240



**Allows surgical instruments to be stored and transported**

## Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019444	HS Code	90189084
Warranty (months)	12	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	24	Bare Unit / Width (cm)	24
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	11.6
Packed unit / Length (cm)	24	Bare Unit / Length (cm)	24
Packed unit / Weight	0.241	Bare Unit / Weight (kg)	0.241

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

 07/04/2026 02:38:13 – For information purposes only, can be subject to change