

**S/S KIDNEY DISH 16 X 8.5 X 3.5CM**

Ref. AI16124



07/04/2026 03:59:15 – For information purposes only, can be subject to change

# S/S KIDNEY DISH 16 X 8.5 X 3.5CM

Ref. [AI16124](#)



**Ideal for storing instruments or temporarily installing medical waste**

## Description

- Decontamination with detergent and disinfectant -
- Autoclavable (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

## Features

|                                    |               |                         |          |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking                         | CE            | Medical device          | Class I  |
| Bar code                           | 3661040019468 | HS Code                 | 90189084 |
| Warranty (months)                  | 12            | Country                 | Pakistan |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                         |          |
| Packed unit / Width (cm)           | 17            | Bare Unit / Width (cm)  | 8.5      |
| Packed unit / Height (cm)          | 7             | Bare Unit / Height (cm) | 3.5      |
| Packed unit / Length (cm)          | 21            | Bare Unit / Length (cm) | 16       |
| Packed unit / Weight               | 0.077         | Bare Unit / Weight (kg) | 0.077    |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                         |          |
| Minimum temperature                | NA            | Max temperature         | NA       |
| Minimum humidity                   | NA            | Maximum humidity        | NA       |
| <b>OTHER FEATURES</b>              |               |                         |          |
| Packed unit / Quantity             | 1             |                         |          |

📅 07/04/2026 03:59:15 – For information purposes only, can be subject to change