

**S/S BOX PERFORATED 25 X 12 X 6 CM
PERFORATED TOP / BOTTOM**

Ref. AI25126P



07/04/2026 02:40:48 – For information purposes only, can be subject to change

S/S BOX PERFORATED 25 X 12 X 6 CM PERFORATED TOP / BOTTOM



Ref. [AI25126P](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067773	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Bare Unit / Width (cm)	12.9
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	6.6
Packed unit / Length (cm)	26	Bare Unit / Length (cm)	25.7
Packed unit / Weight	0.67	Bare Unit / Weight (kg)	0.613
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 02:40:48 – For information purposes only, can be subject to change