

**S/S BOX PERFORATED 30 X 12 X 6 CM
PERFORATED TOP / BOTTOM**

Ref. AI30126P



07/04/2026 01:10:51 – For information purposes only, can be subject to change

**S/S BOX PERFORATED 30 X 12 X 6 CM
PERFORATED TOP / BOTTOM**



Ref. AI30126P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067827	HS Code	90189084
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	32	Bare Unit / Length (cm)	30
Packed unit / Weight	0.67	Bare Unit / Weight (kg)	0.67

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

07/04/2026 01:10:51 – For information purposes only, can be subject to change