

**S/S BOX 30 X 12 X 6 CM**

Ref. AI30126



07/04/2026 00:21:40 – For information purposes only, can be subject to change

# S/S BOX 30 X 12 X 6 CM

Ref. [AI30126](#)



**Convenient for storing surgical instruments**

## Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019598	HS Code	90189084
Warranty (months)	12	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	12	Bare Unit / Width (cm)	12
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	30
Packed unit / Weight	0.81	Bare Unit / Weight (kg)	0.758

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 00:21:40 – For information purposes only, can be subject to change