

# REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 3 MM), 15CM



Ref. IA08409

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040006406	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	18
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.015	Bare Unit / Weight (kg)	0.015
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

07/04/2026 00:39:00 – For information purposes only, can be subject to change