

VAGINAL EXAMINATION SET

Ref. IBO103



06/04/2026 22:59:23 – For information purposes only, can be subject to change

VAGINAL EXAMINATION SET

Ref. IBO103



Description

- Box of composition of instruments of general surgery
- Specialty: Gynecology

This box contains:

- Stainless steel box, 30 x 12 x 6 cm
- Dressing forceps, Longuette, straight, 24 cm -IS, Holtex
- Museux forceps, straight, 24 cm, 7 x 5 mm
- Bladder catheter, metal, female, 15 cm
- Collin vaginal speculum, 30 mm, length 105 mm
- Collin vaginal speculum, 40 mm, length 110 mm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040081601	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.6
Bare Unit / Weight (kg)	1.6		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 22:59:23 – For information purposes only, can be subject to change