

INTRA OCULAR LENS IMPLANTATION SET

Ref. IBOI420



Description

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040082837	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.079
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 07/04/2026 03:14:37 – For information purposes only, can be subject to change