

# LISTER BANDAGE SCISSORS 16 CM

Ref. IC06616



06/04/2026 22:51:32 – For information purposes only, can be subject to change

# LISTER BANDAGE SCISSORS 16 CM

Ref. IC06616



**Ideal for small care, especially to cut bandages and dressings without hurting the patient**

## Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040009063	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.059	Bare Unit / Weight (kg)	0.059
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 06/04/2026 22:51:32 – For information purposes only, can be subject to change