

ADSON HOOK SHARP 20 CM

Ref. IC19201



06/04/2026 22:35:19 – For information purposes only, can be subject to change

ADSON HOOK SHARP 20 CM

Ref. IC19201



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040010588	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20.2
Packed unit / Weight	0.027	Bare Unit / Weight (kg)	0.027

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 06/04/2026 22:35:19 – For information purposes only, can be subject to change