

# DOYEN MOUTH-GAG CHILD, 12 CM

Ref. ID03600



CE



07/04/2026 03:14:58 – For information purposes only, can be subject to change

# DOYEN MOUTH-GAG CHILD, 12 CM

Ref. ID03600



## Description

Doyen children's mouth opener in stainless steel with a total length of 12 cm (Hooks: 45 mm - Ends: 21 mm).

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040013442	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	12
Packed unit / Weight	0.081	Bare Unit / Weight (kg)	0.081
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 07/04/2026 03:14:58 – For information purposes only, can be subject to change