

**PROBE N°1**

Ref. IDSOND01



**Description**

Rhein probes

9 different models are available, contact us

06/04/2026 23:29:40 – For information purposes only, can be subject to change

# PROBE N°1

Ref. IDSOND01



## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040016009	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15.5
Packed unit / Weight	0.027	Bare Unit / Weight (kg)	0.027
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

06/04/2026 23:29:40 – For information purposes only, can be subject to change