

**PROBE N°23**

Ref. IDSOND23



**Description**

Rhein probes

9 different models are available, contact us

06/04/2026 23:22:35 – For information purposes only, can be subject to change

# PROBE N°23

Ref. IDSOND23



## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040016085	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.7
Packed unit / Weight	0.023	Bare Unit / Weight (kg)	0.023

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

06/04/2026 23:22:35 – For information purposes only, can be subject to change