

# LABORDE RETRACTOR 3 BRANCHS, 15CM

Ref. IE24230



## Features

|                                    |               |                         |          |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking                         | CE            | Medical device          | Class I  |
| Bar code                           | 3661040017006 | HS Code                 | 90189084 |
| Warranty (months)                  | 24            | Country                 | Pakistan |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                         |          |
| Packed unit / Width (cm)           | nc            | Bare Unit / Width (cm)  | 8.7      |
| Packed unit / Height (cm)          | nc            | Bare Unit / Height (cm) | 3        |
| Packed unit / Length (cm)          | nc            | Bare Unit / Length (cm) | 14.6     |
| Packed unit / Weight               | 0.062         | Bare Unit / Weight (kg) | 0.062    |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                         |          |
| Minimum temperature                | NA            | Max temperature         | NA       |
| Minimum humidity                   | NA            | Maximum humidity        | NA       |
| <b>OTHER FEATURES</b>              |               |                         |          |
| Packed unit / Quantity             | 1             |                         |          |

07/04/2026 03:28:13 – For information purposes only, can be subject to change