

TUFFIER RETRACTOR 50 X 40 MM

Ref. IE29540



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040017167	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	12.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23.8
Packed unit / Weight	0.465	Bare Unit / Weight (kg)	0.462
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 03:32:58 – For information purposes only, can be subject to change