

# YASARGIL/LEYLA SELF-RETAINING RETRACTOR FOR SKULL MOUNTING

Ref. IE32200



06/04/2026 23:06:37 – For information purposes only, can be subject to change

# YASARGIL/LEYLA SELF-RETAINING RETRACTOR FOR SKULL MOUNTING



Ref. IE32200

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040017303	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.45
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

06/04/2026 23:06:37 – For information purposes only, can be subject to change