

COLLIN MALLET DOUBLE HEADS 20 CM 200 GR (STAINLESS STEEL)



Ref. [IM01220](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040020372	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	4	Bare Unit / Width (cm)	3.4
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.8
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	20.7
Packed unit / Weight	0.026	Bare Unit / Weight (kg)	0.026
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 22:43:33 – For information purposes only, can be subject to change