

SPLINT FORCEPS 8 CM WITH MAGNIFIER (DIAM 30 MM)



Ref. IP01400

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040021362	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.9
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	8
Packed unit / Weight	0.02	Bare Unit / Weight (kg)	0.02

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

06/04/2026 23:48:05 – For information purposes only, can be subject to change