

# KOCHER FORCEPS 1X2 TEETH STRAIGHT 15 CM

Ref. IP02015



06/04/2026 23:47:52 – For information purposes only, can be subject to change

# KOCHER FORCEPS 1X2 TEETH STRAIGHT 15 CM



Ref. IP02015

**Allows you to clamber or fix elements (tubing, etc.). Used mainly on sturdy fabrics.**

## Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021393	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15.3
Packed unit / Weight	0.042	Bare Unit / Weight (kg)	0.042
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 06/04/2026 23:47:52 – For information purposes only, can be subject to change