

ROCHESTER-PEAN FORCEPS CURVED 15 CM

Ref. IP02315




CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL



06/04/2026 23:25:36 – For information purposes only, can be subject to change

ROCHESTER-PEAN FORCEPS CURVED 15 CM



Ref. IP02315

Convenient to stop bleeding without damaging fragile tissues during an operation. Also used to perform care (dressings, etc.).

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040085777	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15.5
Packed unit / Weight	0.043	Bare Unit / Weight (kg)	0.043
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 06/04/2026 23:25:36 – For information purposes only, can be subject to change