

KOCHER FORCEPS STRAIGHT 1X2 TH 14 CM (COLLIN LOCK)



Ref. IP02414

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021690	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14
Packed unit / Weight	0.037	Bare Unit / Weight (kg)	0.037
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 01:24:31 – For information purposes only, can be subject to change