

## CHAPUT FORCEPS 13 CM

Ref. IP11613



CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

STAINLESS  
**INOX**  
STEEL

07/04/2026 01:37:21 – For information purposes only, can be subject to change

# CHAPUT FORCEPS 13 CM

Ref. IP11613



## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040022659	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	13
Packed unit / Weight	0.028	Bare Unit / Weight (kg)	0.028
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

07/04/2026 01:37:21 – For information purposes only, can be subject to change