

DOYEN INTESTINAL CLAMP 23 CM CURVED

Ref. IP14523



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

07/04/2026 02:23:10 – For information purposes only, can be subject to change

DOYEN INTESTINAL CLAMP 23 CM CURVED

Ref. IP14523



For transverse occlusion of a section of intestine while suturing, resection or anastomosis

Description

Holtex is specialist in stainless steel surgical instruments since 1980. Holtex instruments are all made of top quality stainless steel in order to guarantee good resistance to corrosion and durability.

> This reference is made with Stainless Steel AISI 410.

The Holtex surgical instruments are classified as medical devices and comply ISO 13485 standards and the 93/42/EC directive.

EN ISO 13485: 2016, ISO 14971: 2019, EN 1041: 2008, and EN ISO 15223-1

Holtex Instruments are guaranteed for 2 years under normal conditions of use. The warranty does not cover wear parts, shocks and maintenance defects made by the user

Packaging and traceability: Holtex instruments are supplied in individual bags, sealed and labelled with : Logo, and address of the manufacturer, the product reference, the detailed designation in French and English, the lot number ("LOT"), the CE marking, and the barcode.

Cleaning: Please read carefully instructions of cleaning and sterilisation of stainless-steel equipment inserted in each individual bag. Instructions can be also downloaded from the product page on our website.

DOYEN INTESTINAL CLAMP 23 CM CURVED

Ref. IP14523



Features

| | | | |
|-------------------|---------------|----------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040022925 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|------|-------------------------|-------|
| Packed unit / Width (cm) | 28 | Bare Unit / Width (cm) | 9 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 1.2 |
| Packed unit / Length (cm) | 9 | Bare Unit / Length (cm) | 28.5 |
| Packed unit / Weight | 0.07 | Bare Unit / Weight (kg) | 0.067 |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

07/04/2026 02:23:10 – For information purposes only, can be subject to change