

KELLY FORCEPS STRAIGHT 16 CM

Ref. IP16416



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

07/04/2026 02:40:31 – For information purposes only, can be subject to change

KELLY FORCEPS STRAIGHT 16 CM

Ref. [IP16416](#)



Used during surgery (stopping bleeding, grabbing compresses, etc.)

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040023038	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.045	Bare Unit / Weight (kg)	0.045

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 02:40:31 – For information purposes only, can be subject to change