

TRIANGULAR SHAPED FORCEPS 16 CM

Ref. IP20416



CE

DISPOSITIF MEDICAL
CLASSE
I_r
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

06/04/2026 23:31:43 – For information purposes only, can be subject to change

TRIANGULAR SHAPED FORCEPS 16 CM

Ref. IP20416



Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040023502	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	15.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.034	Bare Unit / Weight (kg)	0.034

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

06/04/2026 23:31:43 – For information purposes only, can be subject to change