

## MURPHY-PEAN FORCEPS STRAIGHT 15 CM

Ref. IP21415



CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

STAINLESS  
**INOX**  
STEEL

07/04/2026 02:14:12 – For information purposes only, can be subject to change

# MURPHY-PEAN FORCEPS STRAIGHT 15 CM

Ref. [IP21415](#)



**Used during surgical procedures (stop bleeding, hold compress, etc.)**

## Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040023632	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.046		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

📅 07/04/2026 02:14:12 – For information purposes only, can be subject to change