

DOYEN NEEDLE HOLDER 13 CM

Ref. IP30413



07/04/2026 03:19:10 – For information purposes only, can be subject to change

DOYEN NEEDLE HOLDER 13 CM

Ref. IP30413



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040024813	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	13
Packed unit / Weight	0.038	Bare Unit / Weight (kg)	0.038
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 03:19:10 – For information purposes only, can be subject to change