

OLSEN-HÉGAR NEEDLE HOLDER 14 CM

Ref. IP32414



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

06/04/2026 22:54:19 – For information purposes only, can be subject to change

OLSEN-HÉGAR NEEDLE HOLDER 14 CM

Ref. IP32414



Features

| | | | |
|-------------------|---------------|----------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040024936 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|-------|---------------------------|-------|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 14 |
| Packed unit / Weight | 0.034 | Bare Unit / Weight (kg) | 0.034 |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

06/04/2026 22:54:19 – For information purposes only, can be subject to change