

OLSEN-HÉGAR NEEDLE HOLDER 18 CM

Ref. IP32418



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040024943	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.05
Bare Unit / Weight (kg)	0.05		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 01:54:31 – For information purposes only, can be subject to change