

# FINOCHIETTO NEDDLE HOLDER STRAIGHT 15 CM



Ref. IP39415

## Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3661040025285 | HS Code        | 90189084 |
| Warranty (months) | 24            | Country        | Pakistan |

### DIMENSIONS AND WEIGHT

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | nc    | Bare Unit / Width (cm)  | 30    |
| Packed unit / Height (cm) | nc    | Bare Unit / Height (cm) | 0.6   |
| Packed unit / Length (cm) | nc    | Bare Unit / Length (cm) | 15.2  |
| Packed unit / Weight      | 0.045 | Bare Unit / Weight (kg) | 0.045 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

📅 06/04/2026 23:55:10 – For information purposes only, can be subject to change