

# WHITE TONSIL FORCEPS 23 CM

Ref. IP44223



## Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040025506	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23.2
Packed unit / Weight	0.066	Bare Unit / Weight (kg)	0.066
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 07/04/2026 01:40:40 – For information purposes only, can be subject to change