

MEMBRANE PUNCH 19 CM

Ref. IP44401



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040025551	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Bare Unit / Width (cm)	5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.1
Packed unit / Length (cm)	19	Bare Unit / Length (cm)	19
Packed unit / Weight	0.015	Bare Unit / Weight (kg)	0.015
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 01:59:23 – For information purposes only, can be subject to change