

BERGER UTERINE BIOPSY FORCEPS ROUND CUT 20 CM X 5 MM



Ref. IP44720

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040025650	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	10.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	27
Packed unit / Weight	0.085	Bare Unit / Weight (kg)	0.085
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 22:30:44 – For information purposes only, can be subject to change