

SCHAEDEL TOWEL FORCEPS 9 CM

Ref. IP47009



CE



Schaedel towel forceps 9 cm

Description

Schaedel towel forceps 9 cm

07/04/2026 03:43:58 – For information purposes only, can be subject to change

SCHAEDEL TOWEL FORCEPS 9 CM

Ref. IP47009



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040026749	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.9
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	9
Packed unit / Weight	0.016	Bare Unit / Weight (kg)	0.016
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 03:43:58 – For information purposes only, can be subject to change