

STAINLESS STEEL RULER 30 CM

Ref. IR10030



07/04/2026 00:05:36 – For information purposes only, can be subject to change

STAINLESS STEEL RULER 30 CM

Ref. IR10030



Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040029481	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	3	Bare Unit / Width (cm)	3.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.1
Packed unit / Length (cm)	33	Bare Unit / Length (cm)	33.2
Packed unit / Weight	0.017	Bare Unit / Weight (kg)	0.017

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 00:05:36 – For information purposes only, can be subject to change