

# HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°7/8

Ref. IS01409



07/04/2026 02:06:44 – For information purposes only, can be subject to change

# HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°7/8

Ref. IS01409



## Description

Hégar Double Ended Dilator Stainless Steel

15 Available sizes, contact us

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029580	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.9
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	21
Packed unit / Weight	0.048		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

📅 07/04/2026 02:06:44 – For information purposes only, can be subject to change