

# HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°27/28

Ref. IS01427

Exemple de bougies Hégar double, inox.  
15 tailles disponibles, nous consulter.



06/04/2026 23:42:55 – For information purposes only, can be subject to change

# HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°27/28

Ref. IS01427



## Description

Hégar Double Ended Dilator Stainless Steel

15 Available sizes, contact us

## Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3661040029672 | HS Code        | 90189084 |
| Warranty (months) | 24            | Country        | Pakistan |

### DIMENSIONS AND WEIGHT

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | nc    | Bare Unit / Width (cm)  | 6     |
| Packed unit / Height (cm) | nc    | Bare Unit / Height (cm) | 2.7   |
| Packed unit / Length (cm) | nc    | Bare Unit / Length (cm) | 22    |
| Packed unit / Weight      | 0.183 | Bare Unit / Weight (kg) | 0.183 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

📅 06/04/2026 23:42:55 – For information purposes only, can be subject to change