

# HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°7

Ref. IS01607



Exemple de bougies Hégar simple, inox.  
30 tailles disponibles, nous consulter



06/04/2026 23:38:29 – For information purposes only, can be subject to change

# HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°7



Ref. IS01607

## Description

Hégar Single Ended Dilator Stainless Steel

30 Sizes available, contact us

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029771	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.056	Bare Unit / Weight (kg)	0.056

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 06/04/2026 23:38:29 – For information purposes only, can be subject to change