

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°8

Ref. IS01608



CE



Exemple de bougies Hégar simple, inox.
30 tailles disponibles, nous consulter



07/04/2026 00:49:04 – For information purposes only, can be subject to change

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°8



Ref. IS01608

Description

Hégar Single Ended Dilator Stainless Steel

30 Sizes available, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029788	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.8
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.073	Bare Unit / Weight (kg)	0.073
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/04/2026 00:49:04 – For information purposes only, can be subject to change