

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°13

Ref. IS01613



CE



Exemple de bougies Hégar simple, inox.
30 tailles disponibles, nous consulter



07/04/2026 01:05:44 – For information purposes only, can be subject to change

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°13



Ref. IS01613

Description

Hégar Single Ended Dilator Stainless Steel

30 Sizes available, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029832	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19.2
Packed unit / Weight	0.057	Bare Unit / Weight (kg)	0.057

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 01:05:44 – For information purposes only, can be subject to change