

# HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°26

Ref. IS01626



CE



Exemple de bougies Hégar simple, inox.  
30 tailles disponibles, nous consulter



07/04/2026 01:41:53 – For information purposes only, can be subject to change

# HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°26



Ref. IS01626

## Description

Hégar Single Ended Dilator Stainless Steel

30 Sizes available, contact us

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029962	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.153	Bare Unit / Weight (kg)	0.153

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

07/04/2026 01:41:53 – For information purposes only, can be subject to change