

METAL FEMALE CATHETER 15 CM X 4 MM

Ref. IS03400



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040087573	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.006	Bare Unit / Weight (kg)	0.006
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 00:29:20 – For information purposes only, can be subject to change