

COLLIN VAGINAL SPECULUM 35 MM X 110 MM

Ref. IS04635



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2
2 ANS
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

06/04/2026 23:25:09 – For information purposes only, can be subject to change

COLLIN VAGINAL SPECULUM 35 MM X 110 MM



Ref. IS04635

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040030234	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	3.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	11
Packed unit / Weight	0.256	Bare Unit / Weight (kg)	0.256
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 23:25:09 – For information purposes only, can be subject to change