

COLLIN VAGINAL SPECULUM 38 MM X 110 MM

Ref. IS04638



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2
2 ANS
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

07/04/2026 03:41:08 – For information purposes only, can be subject to change

COLLIN VAGINAL SPECULUM 38 MM X 110 MM



Ref. IS04638

Features

| | | | |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040030241 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 3.8 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 6 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 11 |
| Packed unit / Weight | 0.274 | Bare Unit / Weight (kg) | 0.274 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

07/04/2026 03:41:08 – For information purposes only, can be subject to change