

**STAINLESS STEEL INSTRUMENT TRAY 46 X  
31 X 5 CM**

Ref. AI46316



**Ideal for storing small surgical instruments**

07/04/2026 03:22:50 – For information purposes only, can be subject to change

# STAINLESS STEEL INSTRUMENT TRAY 46 X 31 X 5 CM

Ref. AI46316



## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040066868	HS Code	90189084
Warranty (months)	12	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

07/04/2026 03:22:50 – For information purposes only, can be subject to change