

**KOCHER FORCEPS STRAIGHT 1X2 TH 18
CM (COLLIN LOCK)**

Ref. IP02418



06/04/2026 22:48:34 – For information purposes only, can be subject to change

KOCHER FORCEPS STRAIGHT 1X2 TH 18 CM (COLLIN LOCK)



Ref. IP02418

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040106175	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.037
Bare Unit / Weight (kg)	0.037		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 22:48:34 – For information purposes only, can be subject to change